



PGAM2 mouse mAb

Catalog No	YP-mAb-17239
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PGAM2 PGAMM
Immunogen	Synthesized peptide derived from human N-terminal PGAM2
Specificity	This antibody detects endogenous levels of PGAM2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Mouse,Monoclonal
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BPG dependent PGAM 2;BPG-dependent PGAM 2;EC:5.4.2.11;Muscle-specific phosphoglycerate mutase;PGAM M;PGAM2;PGAM2 PGAMM; Phosphoglycerate mutase 2 (EC 3.1.3.13) (EC 5.4.2.1) (EC 5.4.2.4) (BPG-dependent PGAM 2) (Muscle-specific phosphoglycerate mutase) (Phosphoglycerate mutase isozyme M) (PGAM-M)
Observed Band 28kD	29 kDa
Calculated Molecular Weight	29 kDa
Cell Pathway	
Tissue Specificity	Expressed in the heart and muscle. Not found in the liver and brain.
Function	Interconversion of 3- and 2-phosphoglycerate with 2,3-bisphosphoglycerate as the primer of the reaction. Can also catalyze the reaction of EC 5.4.2.4 (synthase), but with a reduced activity.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images