



G0S2 mouse mAb

Catalog No	YP-mAb-17243
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	G0S2
Immunogen	Synthesized peptide derived from human C-terminal G0S2
Specificity	This antibody detects endogenous levels of G0S2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Mouse,Monoclonal
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GS 2;G/G1 switch protein 2;G0S2;G0/G1 switch protein 2 (G0/G1 switch regulatory protein 2) (Putative lymphocyte G0/G1 switch gene)
Observed Band	11kDa
Calculated Molecular Weight	11 kDa
Cell Pathway	Mitochondrion .
Tissue Specificity	Widely expressed with highest levels in peripheral blood, skeletal muscle and heart, followed by kidney and liver.
Function	Promotes apoptosis by binding to BCL2, hence preventing the formation of protective BCL2-BAX heterodimers.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

