



CHCH1 mouse mAb

Catalog No	YP-mAb-17269
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CHCHD1 C10orf34
Immunogen	Synthesized peptide derived from human N-terminal CHCH1
Specificity	This antibody detects endogenous levels of CHCH1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Mouse, Monoclonal
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CHCHD1;CHCHD1 C10orf34;Coiled-coil-helix-coiled-coil-helix domain-containing protein 1 (Nuclear protein C2360)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion . Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images