



PRD16 mouse mAb

Catalog No	YP-mAb-17279
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PRDM16 KIAA1675 MEL1 PFM13
Immunogen	Synthesized peptide derived from human C-terminal PRD16
Specificity	This antibody detects endogenous levels of PRD16 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Mouse, Monoclonal
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PRDM16;PRDM16 KIAA1675 MEL1 PFM13;PR domain zinc finger protein 16 (PR domain-containing protein 16) (Transcription factor MEL1) (MDS1/EVI1-like gene 1)
Observed Band	140kD
Calculated Molecular Weight	140kD
Cell Pathway	Nucleus . Cytoplasm .
Tissue Specificity	Expressed in uterus and kidney. Expressed in both cardiomyocytes and interstitial cells.
Function	Binds DNA and functions as a transcriptional regulator . Displays histone methyltransferase activity and monomethylates 'Lys-9' of histone H3 (H3K9me1) in vitro (By similarity). Probably catalyzes the monomethylation of free histone H3 in the cytoplasm which is then transported to the nucleus and incorporated into nucleosomes where SUV39H methyltransferases use it as a substrate to catalyze histone H3 'Lys-9' trimethylation (By similarity). Likely to be one of the primary histone methyltransferases along with MECOM/PRDM3 that direct cytoplasmic H3K9me1 methylation (By similarity). Functions in the differentiation of brown adipose tissue (BAT) which is specialized in dissipating chemical energy in the form of heat in response to cold or excess feeding while white



adipose tissue (WAT) is specialized in the storage of excess energy and the control of systemic metabolism (By similarity). To

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images