



ADX mouse mAb

Catalog No	YP-mAb-17285
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FDX1 ADX
Immunogen	Synthesized peptide derived from human N-terminal ADX
Specificity	This antibody detects endogenous levels of ADX at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Mouse,Monoclonal
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ADX;FDX;Ferredoxin-1;FDX1;FDX1 ADX;Adrenodoxin;mitochondrial (Adrenal ferredoxin) (Ferredoxin-1) (Hepatoredoxin); mitochondrial (Adrenal ferredoxin) (Ferredoxin-1) (Hepatoredoxin)
Observed Band	14 kDa
Calculated Molecular Weight	184 aa, 19 kDa
Cell Pathway	Mitochondrion matrix .
Tissue Specificity	Highest levels in the adrenal gland (at protein level). Also detected in kidney and testis (at protein level).
Function	Essential for the synthesis of various steroid hormones . Participates in the reduction of mitochondrial cytochrome P450 for steroidogenesis . Transfers electrons from adrenodoxin reductase to CYP11A1, a cytochrome P450 that catalyzes cholesterol side-chain cleavage . Does not form a ternary complex with adrenodoxin reductase and CYP11A1 but shuttles between the two enzymes to transfer electrons (By similarity).
Background	



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

