



MTTP/MTP mouse mAb

Catalog No	YP-mAb-17686
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MTTP MTP
Immunogen	Synthesized peptide derived from human MTTP/MTP
Specificity	This antibody detects endogenous levels of MTTP/MTP at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, and 0.3090% sodium azide.
Source	Mouse, Monoclonal
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MTTP;MTP; NDGOA;PARK5;PGP95;SPG79;PGP9.5;UCHL-1;Uch-L1;HEL-117; PGP 9.5;HEL-S-53;MTTP MTP;Microsomal triglyceride transfer protein large subunit
Observed Band	98kD
Calculated Molecular Weight	98kD
Cell Pathway	
Tissue Specificity	Liver and small intestine. Also found in ovary, testis and kidney.
Function	Catalyzes the transport of triglyceride, cholesteryl ester, and phospholipid between phospholipid surfaces . Required for the assembly and secretion of plasma lipoproteins that contain apolipoprotein B . May be involved in regulating cholesteryl ester biosynthesis in cells that produce lipoproteins (By similarity).
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images