



PIP5K1A mouse mAb

Catalog No	YP-mAb-17703
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PIP5K1A
Immunogen	Fusion protein of human PIP5K1A
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	68 kDa type I phosphatidylinositol 4-phosphate 5-kinase alpha;EC:2.7.1.68; Phosphatidylinositol 4-phosphate 5-kinase type I alpha;Phosphatidylinositol 4-phosphate 5-kinase type-1 alpha;PIP5K1 alpha;PIP5K1A;phosphatidylinositol-4-phosphate 5-kinase;type I;alpha
Observed Band	63 kDa
Calculated Molecular Weight	56 kDa, 63 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Sodium-independent, high-affinity exchange of anionic amino acids with high specificity for anionic form of cystine and glutamate.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images