



AXUD1 mouse mAb

Catalog No	YP-mAb-17728
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CSRNP1 AXUD1 TAIP3
Immunogen	Synthesized peptide derived from human AXUD1
Specificity	This antibody detects endogenous levels of AXUD1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AXIN1 up regulated 1;AXUD1;CSRNP1;CSRNP1 Ens;FAM130B;TAIP 3;URAX1; CSRNP1 AXUD1 TAIP3;Cysteine/serine-rich nuclear protein 1 (CSRNP-1) (Axin-1 up-regulated gene 1 protein) (Protein URAX1) (TGF-beta-induced apoptosis protein 3) (TAIP-3)
Observed Band	64-70 kDa
Calculated Molecular Weight	64 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitous. Most abundantly expressed in lung, placenta, skeletal muscle, pancreas and leukocyte. Frequently down-regulated in lung, kidney, liver and colon cancers compared with their corresponding normal tissues.
Function	Binds to the consensus sequence 5'-AGAGTG-3' and has transcriptional activator activity (By similarity). May have a tumor-suppressor function. May play a role in apoptosis.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images