



KIR2DL5A mouse mAb

Catalog No	YP-mAb-17731
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KIR2DL5A
Immunogen	Synthetic peptide of human KIR2DL5A
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CD158F;CD158f1;KI2LA_HUMAN;Killer cell immunoglobulin-like receptor 2DL5A;killer cell immunoglobulin-like receptor KIR2DL5A;killer cell immunoglobulin-like receptor;two domains;long cytoplasmic tail;5;5A;killer-cell Ig-like receptor;KIR2DL5;KIR2DL5.1;KIR2DL5.3;KIR2DL5A;KIR2DL5B;killer cell immunoglobulin like receptor;two Ig domains and long cytoplasmic tail 5A
Observed Band	Calculated MW: 41 kDa; Observed MW: 41 kDa
Calculated Molecular Weight	41 kDa; 41kD(Calculated).
Cell Pathway	
Tissue Specificity	
Function	
Background	Receptor on natural killer (NK) cells for HLA-C alleles.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images