



KCNK2 mouse mAb

Catalog No	YP-mAb-17751
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KCNK2
Immunogen	Synthetic peptide of human KCNK2
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KCNK2;TREK1;K2p2.1;KCNK 2;Outward rectifying potassium channel protein TREK-1;potassium two pore domain channel subfamily K member 2;TREK;TPKC1;TREK-1;hTREK-1c;hTREK-1e
Observed Band	40-50 kDa
Calculated Molecular Weight	411 aa, 46 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Ion channel that contributes to passive transmembrane potassium transport.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images