



CD86 mouse mAb

Catalog No	YP-mAb-17790
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CD86
Immunogen	Recombinant protein of human CD86 (NP_787058.4).
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CD86;CD28LG2;T-lymphocyte activation antigen CD86;Activation B7-2 antigen; B70;BU63;CTLA-4 counter-receptor B7.2;FUN-1;CD antigen CD86
Observed Band	Calculated MW: 38 kDa; Observed MW: 74 kDa
Calculated Molecular Weight	Calculated MW: 38 kDa; Observed MW: 74 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Receptor involved in the costimulatory signal essential for T-lymphocyte proliferation and interleukin-2 production, by binding CD28 or CTLA-4. May play a critical role in the early events of T-cell activation and costimulation of naive T-cells, such as deciding between immunity and anergy that is made by T-cells within 24 hours after activation. Isoform 2 interferes with the formation of CD86 clusters, and thus acts as a negative regulator of T-cell activation.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images