



MTX2 mouse mAb

Catalog No	YP-mAb-17793
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MTX2
Immunogen	Synthesized peptide derived from human MTX2
Specificity	This antibody detects endogenous levels of MTX2 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Metaxin-2;metaxin 2;MTX2;Metaxin-2 (Mitochondrial outer membrane import complex protein 2)
Observed Band	25-29 kDa
Calculated Molecular Weight	263 aa, 30 kDa
Cell Pathway	Mitochondrion outer membrane . Mitochondrion .
Tissue Specificity	
Function	Involved in transport of proteins into the mitochondrion.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

