



ACSF2 mouse mAb

Catalog No	YP-mAb-17824
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ACSF2
Immunogen	Fusion protein of human ACSF2
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UNQ493/PRO1009;Medium-chain acyl-CoA ligase ACSF2;mitochondrial; EC:6.2.1.2;ACSF 2;ACSF2;acyl-CoA synthetase family member 2;ACSMW; AVYV493
Observed Band	55-70 kDa
Calculated Molecular Weight	615 aa, 68 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Acyl-CoA synthases catalyze the initial reaction in fatty acid metabolism, by forming a thioester with CoA.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images