



CHGA mouse mAb

Catalog No	YP-mAb-17833
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CHGA
Immunogen	Synthetic peptide of human CHGA
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CGA;CHGA;AL-11;Chromogranin-A;EA-92;beta Granin;Chromogranin A;ChrA;SP-I;Vasostatin I;Vasostatin-2;Vasostatin II;Pancreastatin
Observed Band	Calculated MW: 51 kDa; Observed MW: 51 kDa
Calculated Molecular Weight	51 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Pancreastatin strongly inhibits glucose induced insulin release from the pancreas.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images