



FOXO6 mouse mAb

Catalog No	YP-mAb-17835
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FOXO6
Immunogen	Synthesized peptide derived from human FOXO6
Specificity	This antibody detects endogenous levels of FOXO6 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	forkhead box protein O6;FOXO6
Observed Band	70-80 kDa
Calculated Molecular Weight	51 kDa
Cell Pathway	Cytoplasm . Nucleus . When phosphorylated, translocated from nucleus to cytoplasm. High nuclear localization after stimulation with growth factors (By similarity). .
Tissue Specificity	
Function	Transcriptional activator.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images