



RPS29 mouse mAb

Catalog No	YP-mAb-17836
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPS29
Immunogen	Synthesized peptide derived from human RPS29
Specificity	This antibody detects endogenous levels of RPS29 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	40S ribosomal protein S29;ribosomal protein S29;RPS29
Observed Band	7 kDa
Calculated Molecular Weight	56 aa, 7 kDa
Cell Pathway	Cytoplasm, cytosol . Cytoplasm . Rough endoplasmic reticulum . Detected on cytosolic polysomes (PubMed:25957688). Detected in ribosomes that are associated with the rough endoplasmic reticulum (By similarity). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

