



REG3G mouse mAb

Catalog No	YP-mAb-17841
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	REG3G
Immunogen	Fusion protein of human REG3G
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Pancreatitis-associated protein 1B;Pancreatitis-associated protein IB;PAP IB;PAP1B;PAP-1B;REG3G;regenerating islet-derived 3 gamma;PAPIB;UNQ429;LPPM429;REG III;REG-III
Observed Band	18 kDa
Calculated Molecular Weight	19 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Bactericidal C-type lectin which acts exclusively against Gram-positive bacteria and mediates bacterial killing by binding to surface-exposed carbohydrate moieties of peptidoglycan.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images