



CD62P mouse mAb

Catalog No	YP-mAb-17859
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SELP
Immunogen	The antiserum was produced against synthesized peptide derived from the N-terminal region of human SELP.
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	P-selectin / CD62P;CD62P;P-selectin;AK4;CD62 antigen-like family member P; SELP;GMRP;GRMP;Granule membrane protein 140;GMP-140;Leukocyte-endothelial cell adhesion molecule 3;LECAM3;Platelet activation dependent granule-external membrane protein;PADGEM
Observed Band	Calculated MW: 91 kDa; Observed MW: 91 kDa
Calculated Molecular Weight	Calculated MW: 91 kDa; Observed MW: 91 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Ca(2+)-dependent receptor for myeloid cells that binds to carbohydrates on neutrophils and monocytes. Mediates the interaction of activated endothelial cells or platelets with leukocytes. The ligand recognized is sialyl-Lewis X. Mediates rapid rolling of leukocyte rolling over vascular surfaces during the initial steps in inflammation through interaction with PSGL1.

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images