



CCR6 mouse mAb

Catalog No	YP-mAb-17863
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CCR6
Immunogen	Synthetic peptide of human CCR6
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C C chemokine receptor like 2;C-C chemokine receptor-like 2;CCR11;CCR6; Chemokine receptor CCR11; CCR6;CKRL3;CMKBR6;GPR29;STRL22;C-C chemokine receptor type 6;C-C CKR-6;CC-CKR-6;CCR-6;Chemokine receptor-like 3;CKR-L3;DRY6;G-protein coupled receptor 29;GPR-CY4;GPRCY4;LARC receptor;CD196
Observed Band	Calculated MW: 42 kDa; Observed MW: 42 kDa
Calculated Molecular Weight	Calculated MW: 42 kDa; Observed MW: 42 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Receptor for the C-C type chemokine CCL20.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images