



MOG mouse mAb

Catalog No	YP-mAb-17890
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MOG
Immunogen	A synthesized peptide derived from human Myelin oligodendrocyte glycoprotein
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MOGIG 2;Myelin-oligodendrocyte glycoprotein;MOG;BTN6;BTNL11;MOGIG2;NRCLP7
Observed Band	Calculated MW: 28 kDa; Observed MW: 28 kDa
Calculated Molecular Weight	Calculated MW: 28 kDa; Observed MW: 28 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Mediates homophilic cell-cell adhesion (By similarity). Minor component of the myelin sheath. May be involved in completion and/or maintenance of the myelin sheath and in cell-cell communication.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images