



CREST mouse mAb

Catalog No	YP-mAb-17901
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SS18L1 CREST KIAA0693
Immunogen	Synthesized peptide derived from human CREST
Specificity	This antibody detects endogenous levels of CREST at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SS18L1;Calcium-responsive transactivator;KIAA0693;LP2261;SS18 like protein 1;SS18L1 CREST KIAA0693;Calcium-responsive transactivator (SS18-like protein 1) (SYT homolog 1)
Observed Band	55 kDa
Calculated Molecular Weight	396 aa, 43 kDa
Cell Pathway	Nucleus . Chromosome, centromere, kinetochore . Localizes to nuclear bodies. Colocalizes with SGO1 at kinetochore (By similarity). .
Tissue Specificity	Ubiquitous; with lowest levels in spleen.
Function	Transcriptional activator which is required for calcium-dependent dendritic growth and branching in cortical neurons. Recruits CREB-binding protein (CREBBP) to nuclear bodies. Component of the CREST-BRG1 complex, a multiprotein complex that regulates promoter activation by orchestrating a calcium-dependent release of a repressor complex and a recruitment of an activator complex. In resting neurons, transcription of the c-FOS promoter is inhibited by BRG1-dependent recruitment of a phospho-RB1-HDAC1 repressor complex. Upon calcium influx, RB1 is dephosphorylated by calcineurin, which leads to release of the repressor complex. At the same time, there is increased recruitment of CREBBP to the promoter by a CREST-dependent mechanism,



which leads to transcriptional activation. The CREST-BRG1 complex also binds to the NR2B promoter, and activity-dependent induction of NR2B expression involves

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images