



DEPDC1 mouse mAb

Catalog No	YP-mAb-17905
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DEPDC1
Immunogen	Fusion protein of human DEPDC1
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	583484J8Rik;Cell cycle control protein SDP35;DEP domain containing 1;DEP domain containing protein 1A;DEP domain-containing protein 1A;DEP.8;DEP1A_HUMAN;DEPDC 1;DEPDC 1A;DEPDC1;DEPDC1-V2;Depdc1a;FLJ2354;SDP 35;SDP35
Observed Band	Calculated MW: 93 kDa; Observed MW: 93 kDa
Calculated Molecular Weight	93kD(Calculated).
Cell Pathway	
Tissue Specificity	
Function	
Background	May be involved in transcriptional regulation as a transcriptional corepressor.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images