



GPR89A/GPR89B mouse mAb

Catalog No	YP-mAb-17906
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GPR89A/GPR89B
Immunogen	Synthetic peptide of human GPR89A/GPR89B
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GPR89A/GPR89B;G protein-coupled receptor 89A/B;GPHR;GPR89;SH120;GPR89B;GPR89C;UNQ192
Observed Band	Calculated MW: 53 kDa; Observed MW: 53 kDa
Calculated Molecular Weight	53kD(Calculated).
Cell Pathway	
Tissue Specificity	
Function	
Background	Voltage dependent anion channel required for acidification and functions of the Golgi apparatus that may function in counter-ion conductance. Plays a role in lymphocyte development, probably by acting as a RABL3 effector in hematopoietic cells.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images