



IL-7 mouse mAb

Catalog No	YP-mAb-17925
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	IL7
Immunogen	Synthetic peptide of human IL7
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IL-7;IL7;interleukin 7
Observed Band	Calculated MW: 20 kDa; Observed MW: 20 kDa
Calculated Molecular Weight	Calculated MW: 20 kDa; Observed MW: 20 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	The protein encoded by this gene is a cytokine important for B and T cell development. This cytokine and the hepatocyte growth factor (HGF) form a heterodimer that functions as a pre-pro-B cell growth-stimulating factor. This cytokine is found to be a cofactor for V(D)J rearrangement of the T cell receptor beta (TCRB) during early T cell development. This cytokine can be produced locally by intestinal epithelial and epithelial goblet cells, and may serve as a regulatory factor for intestinal mucosal lymphocytes. Knockout studies in mice suggested that this cytokine plays an essential role in lymphoid cell survival. Alternative splicing results in multiple transcript variants encoding distinct isoforms. Additional splice variants have been described but their presence in



normal tissues has not been confirmed.

**matters needing
attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images