



MIOX mouse mAb

Catalog No	YP-mAb-17944
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MIOX ALDRL6 KSP32 RSOR
Immunogen	Synthesized peptide derived from human MIOX
Specificity	This antibody detects endogenous levels of MIOX at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Aldehyde reductase like 6;ALDRL6;Inositol oxygenase;Kidney specific protein 32;KSP32;MI oxygenase;MIOX;myo inositol oxygenase;Renal specific oxidoreductase;RSOR;MIOX ALDRL6 KSP32 RSOR;Inositol oxygenase (EC 1.13.99.1) (Aldehyde reductase-like 6) (Kidney-specific protein 32) (Myo-inositol oxygenase) (MI oxygenase) (Renal-specific oxidoreductase)
Observed Band	33 kDa
Calculated Molecular Weight	285 aa, 33 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Kidney specific.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images