



TIMM8A mouse mAb

Catalog No	YP-mAb-17945
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TIMM8A DDP DDP1 TIM8A
Immunogen	Synthesized peptide derived from human TIMM8A
Specificity	This antibody detects endogenous levels of TIMM8A at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TIM8A;DDP;DDP1;DFN1;MTS;TIMM8A;TIMM8A DDP DDP1 TIM8A; Mitochondrial import inner membrane translocase subunit Tim8 A (Deafness dystonia protein 1) (X-linked deafness dystonia protein)
Observed Band	11 kDa
Calculated Molecular Weight	11 kDa
Cell Pathway	Mitochondrion inner membrane ; Peripheral membrane protein ; Intermembrane side .
Tissue Specificity	Highly expressed in fetal and adult brain, followed by fetal lung, liver and kidney. Also expressed in heart, placenta, lung, liver, kidney, pancreas, skeletal muscle and heart.
Function	Mitochondrial intermembrane chaperone that participates in the import and insertion of some multi-pass transmembrane proteins into the mitochondrial inner membrane. Also required for the transfer of beta-barrel precursors from the TOM complex to the sorting and assembly machinery (SAM complex) of the outer membrane. Acts as a chaperone-like protein that protects the hydrophobic precursors from aggregation and guide them through the mitochondrial intermembrane space. The TIMM8-TIMM13 complex mediates the import of proteins such as TIMM23, SLC25A12/ARALAR1 and SLC25A13/ARALAR2,



while the predominant TIMM9-TIMM10 70 kDa complex mediates the import of much more proteins. Probably necessary for normal neurologic development.

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images