



RAB3D mouse mAb

Catalog No	YP-mAb-17952
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RAB3D GOV RAB16
Immunogen	Synthesized peptide derived from human RAB3D
Specificity	This antibody detects endogenous levels of RAB3D at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GOV;RAB16;RAD3D;Ras related protein Rab 3D;Ras-related protein Rab-3D; RAB3D;RAB3D GOV RAB16
Observed Band	24-27 kDa
Calculated Molecular Weight	219 aa, 24 kDa
Cell Pathway	Cell membrane ; Lipid-anchor ; Cytoplasmic side .
Tissue Specificity	Highly expressed in granulocytes of peripheral blood. Constitutively expressed at low levels in all hematopoietic cell lines investigated.
Function	Protein transport. Probably involved in regulated exocytosis (By similarity).
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images