



AXIN2 mouse mAb

Catalog No	YP-mAb-17954
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	AXIN2
Immunogen	Synthesized peptide derived from human AXIN2
Specificity	This antibody detects endogenous levels of AXIN2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Axin-2;Axil;axin 2;Axin like protein;Axin-like protein;AXIN2;Axin-2 (Axin-like protein) (Axil) (Axis inhibition protein 2) (Conductin)
Observed Band	94-100 kDa
Calculated Molecular Weight	843 aa, 94 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Expressed in brain and lymphoblast.
Function	Inhibitor of the Wnt signaling pathway. Down-regulates beta-catenin. Probably facilitate the phosphorylation of beta-catenin and APC by GSK3B.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images