



RNF38 mouse mAb

Catalog No	YP-mAb-17962
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RNF38
Immunogen	Synthesized peptide derived from human RNF38
Specificity	This antibody detects endogenous levels of RNF38 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ring finger protein 38;RNF38
Observed Band	43-58 kDa
Calculated Molecular Weight	515 aa, 58 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Widely expressed with highest levels in testis.
Function	Acts as an E3 ubiquitin-protein ligase able to ubiquitinate p53/TP53 which promotes its relocalization to discrete foci associated with PML nuclear bodies. Exhibits preference for UBE2D2 as a E2 enzyme.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images