



KLHL14 mouse mAb

Catalog No	YP-mAb-17967
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KLHL14 KIAA1384
Immunogen	Synthesized peptide derived from human KLHL14
Specificity	This antibody detects endogenous levels of KLHL14 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein interactor of Torsin-1A;Printor;Protein interactor of torsinA;KIAA1384; Kelch-like protein 14;KLHL14;KLHL14 KIAA1384;Kelch-like protein 14 (Protein interactor of Torsin-1A) (Printor) (Protein interactor of torsinA)
Observed Band	35 kDa, 70 kDa
Calculated Molecular Weight	44 kDa, 71 kDa
Cell Pathway	Cytoplasm, cytosol . Endoplasmic reticulum membrane . Colocalizes with TOR1A at the endoplasmic reticulum level.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images