



RPL6 mouse mAb

Catalog No	YP-mAb-17970
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPL6 TXREB1
Immunogen	Synthesized peptide derived from human RPL6
Specificity	This antibody detects endogenous levels of RPL6 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	60S ribosomal protein L6;Large ribosomal subunit protein eL6;Neoplasm related protein C140;Neoplasm-related protein C140;ribosomal protein L6;RPL6;RPL6 TXREB1;60S ribosomal protein L6 (Neoplasm-related protein C140) (Tax-responsive enhancer element-binding protein 107) (TaxREB107)
Observed Band	34 kDa
Calculated Molecular Weight	33 kDa
Cell Pathway	Cytoplasm, cytosol . Cytoplasm . Rough endoplasmic reticulum . Detected on cytosolic polysomes (PubMed:25957688). Detected in ribosomes that are associated with the rough endoplasmic reticulum (By similarity). .
Tissue Specificity	
Function	Component of the large ribosomal subunit. ; (Microbial infection) Specifically binds to domain C of the Tax-responsive enhancer element in the long terminal repeat of HTLV-I .
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images