



KLHL10 mouse mAb

Catalog No	YP-mAb-17989
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KLHL10
Immunogen	Synthesized peptide derived from human KLHL10
Specificity	This antibody detects endogenous levels of KLHL10 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Kelch-like protein 10;Kelch like protein 10;KLHL10
Observed Band	65-70 kDa
Calculated Molecular Weight	608 aa, 69 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	May be a substrate-specific adapter of a CUL3-based E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins during spermatogenesis.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images