



Syntaxin 11 mouse mAb

Catalog No	YP-mAb-18003
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	STX11
Immunogen	Synthesized peptide derived from human Syntaxin 11
Specificity	This antibody detects endogenous levels of Syntaxin 11 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Syntaxin-11;STX11
Observed Band	33-35 kDa
Calculated Molecular Weight	287 aa, 33 kDa
Cell Pathway	Membrane ; Peripheral membrane protein . Golgi apparatus, trans-Golgi network membrane ; Peripheral membrane protein .
Tissue Specificity	
Function	SNARE that acts to regulate protein transport between late endosomes and the trans-Golgi network.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images