



BTN3A2 mouse mAb

Catalog No	YP-mAb-18004
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BTN3A2 BT3.2 BTF3 BTF4
Immunogen	Synthesized peptide derived from human BTN3A2
Specificity	This antibody detects endogenous levels of BTN3A2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BT3.2;BT3.3;BTF3;BTF4;BTN3A2;BTN3A2 BT3.2 BTF3 BTF4;Butyrophilin subfamily 3 member A2
Observed Band	36 kDa
Calculated Molecular Weight	36 kDa
Cell Pathway	Cell membrane ; Single-pass type I membrane protein .
Tissue Specificity	Detected in T-cells and natural killer cells.
Function	Plays a role in T-cell responses in the adaptive immune response. Inhibits the release of IFNG from activated T-cells.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images