



ARMC10 mouse mAb

Catalog No	YP-mAb-18016
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ARMC10 SVH PSEC0198
Immunogen	Synthesized peptide derived from human ARMC10
Specificity	This antibody detects endogenous levels of ARMC10 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	armadillo repeat containing 10;Armadillo repeat-containing protein 10;PNAS 112; PSEC0198;SVH;ARMC10;ARMC10 SVH PSEC0198;Armadillo repeat-containing protein 10 (Splicing variant involved in hepatocarcinogenesis protein)
Observed Band	31-37 kDa
Calculated Molecular Weight	343 aa, 38 kDa
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass membrane protein .
Tissue Specificity	Expressed in all tissues tested with higher expression in placenta, liver, kidney, heart and brain.
Function	May play a role in cell survival and cell growth. May suppress the transcriptional activity of p53/TP53.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images