



## SSX4 mouse mAb

|                             |                                                                                                                                   |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Catalog No                  | YP-mAb-18020                                                                                                                      |
| Isotype                     | IgG                                                                                                                               |
| Reactivity                  | Human;Mouse;Rat                                                                                                                   |
| Applications                | WB                                                                                                                                |
| Gene Name                   | SSX4 SSX4A; SSX4B                                                                                                                 |
| Immunogen                   | Synthesized peptide derived from human SSX4                                                                                       |
| Specificity                 | This antibody detects endogenous levels of SSX4 at Human                                                                          |
| Formulation                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                           |
| Source                      | Monoclonal, Mouse,IgG                                                                                                             |
| Purification                | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.              |
| Dilution                    | WB 1:500-1:2000                                                                                                                   |
| Concentration               | 1 mg/ml                                                                                                                           |
| Purity                      | ≥90%                                                                                                                              |
| Storage Stability           | -20°C/1 year                                                                                                                      |
| Synonyms                    | Cancer/testis antigen 5.4;CT5.4;Protein SSX4;SSX4;SSX4A;SSX4B;SSX4 SSX4A;Protein SSX4 (Cancer/testis antigen 5.4) (CT5.4)         |
| Observed Band               | 22 kDa                                                                                                                            |
| Calculated Molecular Weight | 188 aa, 22 kDa                                                                                                                    |
| Cell Pathway                |                                                                                                                                   |
| Tissue Specificity          |                                                                                                                                   |
| Function                    | Could act as a modulator of transcription.                                                                                        |
| Background                  |                                                                                                                                   |
| matters needing attention   | Avoid repeated freezing and thawing!                                                                                              |
| Usage suggestions           | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



## Products Images