



SPANXB1 mouse mAb

Catalog No	YP-mAb-18021
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SPANXB1 SPANXB; SPANXB2; SPANXF1
Immunogen	Synthesized peptide derived from human SPANXB1
Specificity	This antibody detects endogenous levels of SPANXB1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	B1;Cancer/testis antigen 11.2;CT11.2;Nuclear-associated protein SPAN-Xb; SPANX B;SPANXB1;SPANXB1 SPANXB;SPANXB2;SPANXF1;Sperm protein associated with the nucleus on the X chromosome B/F (Cancer/testis antigen 11.2) (CT11.2) (Nuclear-associated protein SPAN-Xb) (SPANX-B) (Nuclear-associated protein SPAN-Xf) (SPANX-F)
Observed Band	22 kDa
Calculated Molecular Weight	103 aa, 12 kDa
Cell Pathway	Cytoplasm . Nucleus . Associated with nuclear craters.
Tissue Specificity	Detected in testis and sperm.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images