



DDX43 mouse mAb

Catalog No	YP-mAb-18022
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DDX43 HAGE
Immunogen	Synthesized peptide derived from human DDX43
Specificity	This antibody detects endogenous levels of DDX43 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cancer/testis antigen 13;CT13;DDX43;DEAD box protein 43;DEAD box protein HAGE;DKFZp434H2114;HAGE;Helical antigen;DDX43 HAGE;Probable ATP-dependent RNA helicase DDX43 (EC 3.6.4.13) (Cancer/testis antigen 13) (CT13) (DEAD box protein 43) (DEAD box protein HAGE) (Helical antigen)
Observed Band	68-73 kDa
Calculated Molecular Weight	648 aa, 73 kDa
Cell Pathway	
Tissue Specificity	Expressed in testis. Expressed in many tumors of various histological types at a level that is 100-fold higher than the level observed in normal tissues except testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images