



DNAJB8 mouse mAb

Catalog No	YP-mAb-18027
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DNAJB8
Immunogen	Synthesized peptide derived from human DNAJB8
Specificity	This antibody detects endogenous levels of DNAJB8 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DJ6;DNAJB8;DnaJ homolog subfamily B member 8
Observed Band	26 kDa
Calculated Molecular Weight	232 aa, 26 kDa
Cell Pathway	
Tissue Specificity	
Function	Efficient suppressor of aggregation and toxicity of disease-associated polyglutamine proteins.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images