



NXF2 mouse mAb

Catalog No	YP-mAb-18028
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NXF2 TAPL2; NXF2B
Immunogen	Synthesized peptide derived from human NXF2
Specificity	This antibody detects endogenous levels of NXF2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cancer/testis antigen 39;CT39;nuclear RNA export factor 2;NXF2;NXF2B;TAP like protein 2;TAPL 2;TAPL2;NXF2 TAPL2;Nuclear RNA export factor 2 (Cancer/testis antigen 39) (CT39) (TAP-like protein 2) (TAPL-2)
Observed Band	72 kDa
Calculated Molecular Weight	626 aa, 72 kDa
Cell Pathway	Nucleus, nucleoplasm. Cytoplasm. Localized in the nucleoplasm and at the nuclear envelope. Shuttles between the nucleus and the cytoplasm.
Tissue Specificity	Expressed almost exclusively in testis. Also expressed in several cancers.
Function	Involved in the export of mRNA from the nucleus to the cytoplasm.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images