



CCDC83 mouse mAb

Catalog No	YP-mAb-18030
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CCDC83 HSD9
Immunogen	Synthesized peptide derived from human CCDC83
Specificity	This antibody detects endogenous levels of CCDC83 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CCDC83;CCDC83 HSD9;Coiled-coil domain-containing protein 83
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images