



GSTA2 mouse mAb

Catalog No	YP-mAb-18032
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GSTA2 GST2
Immunogen	Synthesized peptide derived from human GSTA2
Specificity	This antibody detects endogenous levels of GSTA2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:2.5.1.18;Glutathione S-transferase A2;GST class-alpha member 2;GST gamma;GST HA subunit 2;GSTA2;GSTA2 GST2;Glutathione S-transferase A2 (EC 2.5.1.18) (GST HA subunit 2) (GST class-alpha member 2) (GST-gamma) (GSTA2-2) (GTH2)
Observed Band	26 kDa
Calculated Molecular Weight	26 kDa
Cell Pathway	Cytoplasm.
Tissue Specificity	Liver.
Function	Catalyzes the conjugation of glutathione to a large variety of electrophilic compounds.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images