



ZNF320 mouse mAb

Catalog No	YP-mAb-18038
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF320
Immunogen	Synthesized peptide derived from human ZNF320
Specificity	This antibody detects endogenous levels of ZNF320 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZFPL;zinc finger protein 320;ZNF320
Observed Band	40-45 kDa
Calculated Molecular Weight	509 aa, 59 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Expressed at high levels in placenta, pancreas, kidney and heart, and at lower levels in liver, skeletal muscle and brain.
Function	May be involved in transcriptional regulation.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images