



SPATC1 mouse mAb

Catalog No	YP-mAb-18057
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SPATC1 SPATA15 SPRN
Immunogen	Synthesized peptide derived from human SPATC1
Specificity	This antibody detects endogenous levels of SPATC1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SPATA15;SPATC1;SPERIOLIN;SPRN;SPATC1 SPATA15 SPRN;Speriolin (Spermatogenesis and centriole-associated protein 1) (Spermatogenesis-associated protein 15) (Spermatogenic cell-specific Cdc20-binding protein)
Observed Band	60 kDa
Calculated Molecular Weight	500 aa, 52 kDa
Cell Pathway	Cytoplasm . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Colocalizes with the centrosomal pericentrin protein PCNT1 (By similarity). Located in the connecting piece of sperm. .
Tissue Specificity	Detected only in testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images