



GPLOW mouse mAb

Catalog No	YP-mAb-18063
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GPLOW GPATC5 GPATCH5 T54
Immunogen	Synthesized peptide derived from human GPLOW
Specificity	This antibody detects endogenous levels of GPLOW at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	T54;Protein T54;GPATCH5;G-patch domain-containing protein 5;G-patch domain and KOW motifs-containing protein;GPLOW;GPLOW GPATC5 GPATCH5 T54;G patch domain and KOW motifs-containing protein (G patch domain-containing protein 5) (Protein MOS2 homolog) (Protein T54)
Observed Band	52 kDa
Calculated Molecular Weight	52 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	RNA-binding protein involved in pre-mRNA splicing.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images