



FADS1 mouse mAb

Catalog No	YP-mAb-18099
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FADS1 FADSD5
Immunogen	Synthesized peptide derived from human FADS1
Specificity	This antibody detects endogenous levels of FADS1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	D5D;TU12;FADS6;FADSD5;LLCDL1;FADS1;FADS1 FADSD5;Fatty acid desaturase 1 (EC 1.14.19.-) (Delta(5) fatty acid desaturase) (D5D) (Delta(5) desaturase) (Delta-5 desaturase)
Observed Band	52kDa
Calculated Molecular Weight	52kDa
Cell Pathway	[Isoform 1]: Endoplasmic reticulum membrane ; Multi-pass membrane protein . Mitochondrion . ; [Isoform 2]: Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Widely expressed, with highest levels in liver, brain, adrenal gland and heart. Highly expressed in fetal liver and brain.
Function	[Isoform 1]: Acts as a front-end fatty acyl-coenzyme A (CoA) desaturase that introduces a cis double bond at carbon 5 located between a preexisting double bond and the carboxyl end of the fatty acyl chain. Involved in biosynthesis of highly unsaturated fatty acids (HUFA) from the essential polyunsaturated fatty acids (PUFA) linoleic acid (LA) (18:2n-6) and alpha-linolenic acid (ALA) (18:3n-3) precursors. Specifically, desaturates dihomo-gamma-linoleoate (DGLA) (20:3n-6) and eicosatetraenoate (ETA) (20:4n-3) to generate arachidonate (AA) (20:4n-6) and eicosapentaenoate (EPA) (20:5n-3), respectively . As a rate limiting enzyme



for DGLA (20:3n-6) and AA (20:4n-6)-derived eicosanoid biosynthesis, controls the metabolism of inflammatory lipids like prostaglandin E2, critical for efficient acute inflammatory response and maintenance of epithelium homeostasis. Contributes to membrane phospholi

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images