



FIBIN mouse mAb

Catalog No	YP-mAb-18109
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FIBIN PSEC0235
Immunogen	Synthesized peptide derived from human FIBIN
Specificity	This antibody detects endogenous levels of FIBIN at Human, Mouse,Ra
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FIBIN;FIBIN PSEC0235;Fin bud initiation factor homolog
Observed Band	23kD
Calculated Molecular Weight	211 aa, 24 kDa
Cell Pathway	Secreted . Golgi apparatus . Endoplasmic reticulum .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images