



FUBP1 mouse mAb

Catalog No	YP-mAb-18110
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FUBP1
Immunogen	Synthesized peptide derived from human FUBP1
Specificity	This antibody detects endogenous levels of FUBP1 at Human,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Far upstream element-binding protein 1;FUSE-binding protein 1;FUBP1;Far upstream element-binding protein 1 (FBP) (FUSE-binding protein 1) (DNA helicase V) (hDH V)
Observed Band	69 kDa
Calculated Molecular Weight	653 aa, 69 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	Regulates MYC expression by binding to a single-stranded far-upstream element (FUSE) upstream of the MYC promoter. May act both as activator and repressor of transcription.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images