



MDH2 mouse mAb

Catalog No	YP-mAb-18120
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MDH2
Immunogen	Synthesized peptide derived from human MDH2
Specificity	This antibody detects endogenous levels of MDH2 at Human, Mouse,Ra
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:1.1.1.37;MDH;MDH 2;MDH2;Malate dehydrogenase;mitochondrial (EC 1.1.1.37); mitochondrial (EC 1.1.1.37)
Observed Band	33 kDa
Calculated Molecular Weight	36 kDa
Cell Pathway	Mitochondrion matrix .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images