



DUSP5 mouse mAb

Catalog No	YP-mAb-18126
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DUSP5 VH3
Immunogen	Synthesized peptide derived from human DUSP5
Specificity	This antibody detects endogenous levels of DUSP5 at Human,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	VH3;HVVH3;EC:3.1.3.16;DUSP;Dual specificity phosphatase 5;DUSP5;DUSP5 VH3;Dual specificity protein phosphatase 5 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity protein phosphatase hVH3)
Observed Band	42 kDa
Calculated Molecular Weight	42kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	Dual specificity protein phosphatase; active with phosphotyrosine, phosphoserine and phosphothreonine residues. The highest relative activity is toward ERK1.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images